Suck Detection Manual

The Suck Detection Program is designed to perform robust and detailed analysis of a data file generated by the Neonur system while being easy to use. Although you may save the data that is generated in the form of an Excel file, the program displays all relevant information completely within the interface of the program.

The Interface

The program's interface consists of three tabs:

Infant Info: Contains fields for general information about the infant as well as a text box for any notes that the user wishes to document. Information recorded in this tab is not used in the program's analysis, but is recorded in the output Excel files.

File	Edit View Insert Tools Desktop Window I	Help	
10	🛃 🎍 🗞 🔍 🤍 🏐 🐙 🖌 - 🗔 🛙		
ë 🖌	<u> </u>		
Infant	Info Params Data		
		Notes :	
	General Information		*
	Gender : MF		
	Age : 0 Weeks		
	Feed Type : Breast/Bottle		
	Gestational Age : 0 Weeks		÷

Params: Contains a Parameters box for the sample rate, correlation coefficient, and cutoff pressure, as well as the "Run Analysis", "Choose Prototype", and "Add Plot to Excel" buttons. This tab also has a graph for the signal and a single row for summary data after analysis.

Paranis Data														
Paramete	rs			1										
Samp	le Rate :	250	Hz	0.9										
Co	rrelation —			0.8										
Coefficient : 0.8			0.7 -											
	Cutoff :	20	mmHg	0.6										
				0.5 -										
	Run A	nalysis		0.4 -										
	Choose I	Prototype		0.3 -										
				0.2										
	Add Plot	to Excel		0.1 -										
				0										
				0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9	
size of signal(s)	number of sucks	mean adjusted max pressure	mean real max pressure	number of bursts	mean burst duration(s)	mean pau duration	ise (s)	mean number of sucks per	mean inter sue interval(s)	:k				
0	0	0	0		0	0	0		0	0				

Data: Displays the full analysis for each suck and burst detected after the program is executed.



Using the Program

To start, click the open folder icon in the top left corner of the toolbar and select the desired data file. This data file may be anywhere on the user's computer. Once the file has been selected, the empty graph in the Params tab will display the signal of the data chosen.



At this point, the user may wish to to change the preset parameters if desired. Additionally, the user may wish to choose a different prototype other than the default if necessary. To run the analysis, press the "Run Analysis" button. The button will then change to "Running..." Once the analysis is complete, the graph will change to reflect the signal processing, the calculated baseline, and the detected sucks in the form of dotted red lines. The data row at the bottom will also fill with the summary data.



Su	uck#	T.O.O	Duration	Adjusted Pmax	Real Pmax	T of Pmax	I.S.I	A.U.C
	1	27.9280	0.5880	3.5406e+04	3.2482e+04	28.2080	0.9240	2.6742e+06
	2	28.5160	0.9800	1.5997e+04	6.9301e+03	29.1320	0.9320	-8.6119e+03
	3	29.4960	1.1680	4.1565e+04	3.2499e+04	30.0640	11	5.0421e+06
	4	40.8040	0.5040	7.7891e+03	1.7130e+04	41.0640	0	1.6503e+06
	5	43.9640	0.8160	5.9090e+03	1.8204e+04	44.3400	-1	3.0820e+06
	Burs	#	# of 5	Sucks	Durat	ion	I.B.I	
	Burs	* # 1	# of \$	Sucks	Durat	ion 2 7360	I.B.I	10 1400
	Burs	: # 1 2	# of \$	5ucks	Durat	ion 2.7360 0.5040	I.B.I	10.1400
	Burs	# 1 2 3	# of !	Sucks 3 1 1	Durat	2.7360 0.5040 0.8160	I.B.I	10.1400 2.6560 -1

The Data tab will also fill with the full analysis of each individual suck and burst detected.

Saving the Results

At this point, the user may wish to save this analysis in the form of an Excel file. To do so, simply click the save icon in the interface's toolbar. The Excel file saved will have two tabs, the first for the data of each individual detected suck and the second for the summary data.

x∎	. 5	~ (÷ - ÷								test - Exc	el							? 📧	- e	×
FILE	E F	IOME	INS	SERT PA	AGE LAYOUT	FORML	LAS DA	ATA RE	VIEW V	EW AI	DD-INS	JMP TEA	M						Brandon L [Pellucky *	
	} ⊪.	Cali	bri	- 11	- A A	==	<i>≫</i> -	🚏 Wrap	Text	Genera	ı	-	i ≠	.	•	>	ΣA	utoSum ×	Ż▼ ₩		
Paste		в	ΙU	• 🛄 •	👌 - 🛕 -	- E B B	i in	📑 Merg	e & Center	- \$ -	% * 5	8 →8 Cond	litional Fo	rmat as Cell	Insert	Delete Form	at 🥒 C	lear •	Sort & Find &	8	
Clipb	oard 5			Font		-	Alian	ment		5	Number	- Forma	styl	able + Styles +		Cells		Edit	ing		~
					C																
A1		٣	1 2	$\langle \checkmark$	Ĵx Suc	k #															~
	Α		В	С	D	E	F	G	н	I	J	К	L	М	N	0	Р	Q	R	s	
1 S	uck #	T.O	.0	Duration	Adjusted I	Real Pmay 1	of Pmax I	.S.I	A.U.C												
2	1	1	27.82	0.86	41548.21	32481.82	28.208	0.924	2343636												
3		2	28.74	0.56	11271.64	6930.068	29.132	0.932	536543.6												
4		3	29.48	1.04	41565.05	32498.61	30.064	8.072	5332260												
5		4	37.84	0.74	1611.541	8535.15	38.136	1.588	1411940												
6	1	5	39.48	0.56	4834.623	13101.18	39.724	1.34	1516423												
7		5	40.76	0.68	7789.115	17130.04	41.064	3.276	2092369												
8		7	43.96	0.86	5908.984	18204.4	44.34	4.536	3219505												
9		8	48.36	1.04	2954.492	22501.84	48.876	-1	5596940												
10										<u>/</u>											
11										Quick A	alysis (Ctrl	+0)									
12										Use the (Juick Analy	sis tool to									
13										quickly a	nd easily ar	nalyze your dat	ta								
14										with son	ne of Excel's	most useful	-								
15										and form	en as charts iulas.	, color-coaing	-								
16													_								
10																					
10																					
20																					
21																					
22																					
23																					
24																					
25																					
4	F		Sheet1	Sheet2	Sheet3	(+)							: [4							
READ	v			-	1 -11010	1 0			_	_		AVERAGE: 34	48736.5824	COUNT: 72	SUM: 223	319141.27	III II	a 🗖 .		-+	100%
•		2				X											(o)	<u>e</u> 🖻 🔋	N P 12 ()	11:32 / 6/17/2	AM 016

x	B 5	~ ¢	÷							test - Excel								? 📧	- 8	×	
F	ILE H	OME	INSERT PA	AGE LAYOUT	FORM	AULAS D	ATA RE	VIEW V	IEW AD	D-INS J	IMP	TEAM						Brandon L	Dellucky *		
	ک ا	Calibri	- 11	- A A	· = =	≡ %·•	🔐 Wrap	Text	Genera		Ŧ	₽	.	•=	>		utoSum 👻	Azy H	b		
Pa	ste 💉	ΒI	<u>u</u> - 🖽 -	<u>⊳</u> - <u>A</u> -	· = =	≡ € #	E 🔛 Merg	e & Center	- \$ - 9	% * 50	.00 •.0	Conditional Fo Formatting * T	rmatas Cell able≖ Styles	Insert	Delete Form	at 🦉 🤇	lear *	Sort & Find Filter * Select	B/		
Clip	board 🕞		Font	1	5	Alig	nment		6 N	ra Number ra Sty			es Cells				Edit		^		
A7 👻 : 🗙			\times \checkmark	∫x Not	tes													~			
	Α	В	С	D	E	F	G	н	I	J	К	L	М	Ν	0	Р	Q	R	S	A	
1	Gender	Age	Feed Type	Gestationa	al Age																
2	M/F		0 Breast/Bo	0																	
3																					
4	size of si	g numbe	er oʻmean adju	mean real	number o	mean bur	mean pau	mean nun	mean inte	r_suck inte	erval(s										
5	49.54	ł	8 1174.837	1513.831	4	2.386667	4.46	2.333333	2.952571												
6	Natas	-																			
/	Notes	-																			
0																					
10																					
11																					
12																					
13																					
14																					
15																					
16																					
17																					
18																					
19																					
20																					
21																					
22																					
23																					
25																					
20	4	She	et1 Sheet2	Sheet3	4)							4								
DEA	DV.	Jone	Jieeu	Sheets	1 (5	,	_	_			_		۱			ATT IS	a m a .		-	100%	
REA																			11.22 A	100%	
C	9 (X											ţ0	🙆 📼 🎽	• 🖻 🖷 🚸	6/17/20	116	

The user may also add the analyzed data plot to a third sheet by clicking the "Add Plot to Excel" button in the Params tab. After clicking, the plot will appear in this manner:



Additionally, one may wish to zoom in on a piece of the graph and save the zoomed graph in the Excel file. To zoom, click the zoom button on the toolbar and select the region you wish to zoom.



To save this new graph to your existing Excel file, click the "Add Plot to Excel" button. The program will prompt you to overwrite your existing file. Confirm the overwrite and the new graph will appear in Sheet 4. You can save as many altered graphs as you like with each one saving in a new sheet.